

**Ethyl Sulfate (25ng/mL) on Synergi Hydro-RP (m/z 124.9 > 80.1)****Column:** Synergi<sup>TM</sup> 2.5  $\mu$ m Hydro-RP 100 Å, LC Column 100 x 3 mm, Ea**Dimensions:** 100 x 3 mm ID**Order No:** 00D-4387-Y0**Elution Type:** Gradient**Eluent A:** 10mM Ammonium formate**Eluent B:** 50:50 Methanol: Acetonitrile

| Gradient Profile: | Step No. | Time (min) | Pct A | Pct B |
|-------------------|----------|------------|-------|-------|
|                   | 1        | 0          | 100   | 0     |
|                   | 2        | 4          | 20    | 80    |
|                   | 3        | 4.1        | 100   | 0     |
|                   | 4        | 7          | 100   | 0     |

**Flow Rate:** 0.5 mL/min**Col. Temp.:** ambient**Detection:** Mass Spectrometer (MS) @ amu (ambient)**Detector Info:** <a target="\_blank" href="https://sciex.com/products/mass-spectrometers?utm\_campaign=2019%20application%20search&utm\_source=phenomenex&utm\_medium=referral">SCIEX</a>

19780

**ANALYTES:****1** EtS (Ethyl sulfate)