

**Miconazole on Lux 5µm Cellulose-3 in PO****Column:** Lux® 5 µm Cellulose-3, LC Column 250 x 4.6 mm, Ea**Dimensions:** 250 x 4.6 mm ID**Order No:** 00G-4493-E0**Elution Type:** Isocratic**Eluent A:** MeOH / IPA (90:10) DEA (0.1%)

| Gradient | Step No. | Time (min) | Pct A |
|----------|----------|------------|-------|
| Profile: | 1        | 0          | 100   |

**Flow Rate:** 1 mL/min**Col. Temp.:** ambient**Detection:** UV-Vis Abs.-Variable Wave.(UV) @ 220 nm (ambient)**Analyst Note:** SecurityGuard™ Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, Lux Cellulose-3 4 x 3.0 mm, 10/Pk Part No.: AJ0-8622

- Holder Part No.: KJ0-4282

**Products used in this application:**

mAU

400

1

**ANALYTES:****1** Miconazole

Retention Time: 5.82 min

**2** Miconazole

Retention Time: 7.58 min

