


Cold-Pressed Orange Oil on ZB1ms 60x0.25x0.25**Column:** Zebron™ ZB-1MS, GC Cap. Column 60 m x 0.25 mm x 0.25 µm, Ea**Phase:** 100% Dimethylpolysiloxane**Dimensions:** 60 meters x 0.25 mm x 0.25 µm**Order No:** 7KG-G011-11**Oven Profile:** 75°C for 4 min to 250°C at 4°C/min for 5 min**Carrier Gas:** Constant Flow Helium, 0.8 mL/min**Injection:** Split 75:1 0.2 µL @ 275°C**Detection:** Refractive Index (0°C)

15580

**Products used in this application:**

Cold-Pressed Orange Oil on ZB1ms 60x0.25x0.25

ANALYTES:

- 1 α -Pinene
- 2 Sabinene
- 3 beta-Myrcene
- 4 Octanal
- 5 3-Carene
- 6 Limonene
- 7 Nonanal
- 8 Linalool
- 9 cis-Limonene oxide
- 10 trans-Limonene oxide
- 11 Citronellal
- 12 α -Terpineol
- 13 Decanal
- 14 Neral
- 15 Carvone
- 16 Geranial
- 17 α -Cubebene
- 18 Dodecanal
- 19 β -Cubenene
- 20 Valencene
- 21 Cadinine
- 22 Nootkatone

