


Propoxyphene by MSD on a ZB-Drug1**Column:** Zebtron[™] ZB-Drug-1, GC Cap. Column 10 m x 0.18 mm x 0.18 µm, Ea**Phase:****Dimensions:** 10 meters x 0.18 mm x 0.18 µm**Order No:** 7CD-G023-08**Oven Profile:** 200°C to 320°C @ 25 °C/min**Carrier Gas:** Constant Flow Helium, 55 cm/sec**Injection:** Split 10:1 1 µL @ 260°C**Detection:** Mass Selective (MSD) (320°C)**Analyst Note:** Analytes are 24ppm
Sim ions:
(Quant ion*) D11-Propoxyphene (*64,198,213): Propoxyphene (*58,193,208)**Products used in this application:**

17803

**ANALYTES:**

- 1 Propoxyphene-d11
- 2 Propoxyphene

