


Test Conditions for Zebron ZB-1plus

Column: Zebron ZB-1plus Capillary GC, Column 30m x 0.25mm x 0.25µm, Ea
Phase: Proprietary Pesticides Phase
Dimensions: 30 meters x 0.25 mm x 0.25 µm
Order No: 7HG-G031-11
Oven Profile: 140°C Isothermal

Carrier Gas: Constant Flow Hydrogen, 1.18 mL/min
Injection: Split 100:1 1 µL @ 250°C
Detection: Flame Ionization (FID) (325°C)
Detector Info: 20mL/min Nitrogen, 30mL/min Hydrogen, 350mL/min Air (total 400mL/min)

Analyst Note: AG0-7805 Zebron ZB-1ms Column Performance Test Mix (ZB-1ms 30m x 0.25mm x 0.25µm). 1.0µl of 500µg/mL

15545

**Products used in this application:****ANALYTES:**

- 1 Decane
- 2 2-Ethylhexanoic acid
- 3 4-Chlorophenol
- 4 Naphthalene
- 5 Tridecane
- 6 1-Undecanol
- 7 Dicyclohexylamine
- 8 Pentadecane

